

ANALOG TENS MACHINE BY METRON: PROTENS

~~\$89.99~~ \$69.99



A regular TENS with adjustable pulse width, pulse rate, intensity, modes and an accurate timer of 15, 30 and 60 minutes, will stop output at preset time. The unit comes with the TENS stimulator unit, carrying case, 2X TENS leads, 9V battery, operating manual and packet of 4 Electrodes. Includes a belt clip.



SKU: S2646

Categories: [Electrotherapy](#), [Pain Management](#), [TENS Machines](#)

PRODUCT DESCRIPTION

The ProTENS unit has a fully adjustable parameter and provides relief from acute or chronic pain such as back pain, sports injury, general muscle pain, etc. Battery powered electrical unit which uses electrodes to deliver electrical pulses to the nerve fibres relieving pain.

ADDITIONAL INFORMATION

Channel	Dual, isolated between channels
Modes	Continuous, Burst, Modulation
Pulse Intensity	Adjustable 0-80mA peak into 500 ohm load each channel, constant current
Pulse Rate	2Hz-150Hz (adjustable)
Pulse Width	30µs-260µs (adjustable)
Time	Continuous, 30 min., 60 min.
Burst Mode	Burst consists 2 burst per sec at 100 Hz
Wave Form	Asymmetrical Bi-Phasic square pulse
Dimensions	95(H) x 60(W) x 23(D) mm
Weight	115 grams (battery included)
Warranty	3 Years
Brand	Metron