



SKU: S7531

KCARE TRANSIT SLIDING FOOTPLATE SHOWER COMMODORE

From: \$1,299.99

The KCare Transit Sliding Footplate Shower Commode offers secure foot support with a sliding design for comfortable, assisted use during showering and toileting.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The KCare Transit Sliding Footplate Shower Commode is a robust and versatile commode to support carer-assisted showering and toileting. It features a sliding footplate that helps to stabilise feet during transfers.

Flip-up armrests assist with side transfers, while the 3-strap backrest enhances comfort.

The steel frame is durable and rust-resistant, suitable for daily wet-area use.

The unit features an integrated push handle and lockable castors to provide safe carer-assisted mobility.

Note: The seat, bed pan and carrier are sold separately.

Who it's for

- For clients seeking a standard-size transit showering commode with extra foot support.
- For therapists recommending a reliable, versatile commode to support daily hygiene tasks.

Key features

- **Transit design** – designed for carer-assisted mobility
- **44.5 cm seat width** – suits a range of individuals while remaining compact

- **Sliding footplate** - offers stable foot support
- **3-strap backrest** - enhances comfort and accessibility
- **Integrated push handle** - allows for easy carer mobility
 - **Flip-up armrests** - support safe, lateral transfers
- **Stainless-steel frame** - durable and corrosion-resistant
- **Four lockable castors** - ensure safe, stable movement

ADDITIONAL INFORMATION

Safe Working Load (SWL)	200 kg
Gross Weight	10 kg
Overall Width	48 cm
Overall Depth	84 cm
Overall Height	100.5 cm
Seat Width	44.5 cm
Seat Depth	43 cm
Seat Height	55 cm
Backrest Height	40 cm
Brand	KCare Healthcare

[KCare Shower Commode User Manual](#)

[KCare Shower Commode Cleaning Instructions](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

No FAQ Found