




**SKU:** S14977

## UNIVERSAL MOBILITY SCOOTER CANOPY

**\$199.99**

The ILS mobility scooter canopy combo is an essential accessory for anyone enjoying some outside time on their scooters. Fits easily to most scooters with a 1 x 1 seat bracket under the rear of their scooter seat. Manufactured high quality water resistant materials and will suit most medium and large scooters.

- Cover adopts high quality Water-resistant Material
- User-friendly collapsible design
- Lightweight metal frame
- Velcro strips attach canopy to hardware
- Adjustable in height and easily removable
- Connecting hardware is included



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

It connects universally to the scooter seat back bracket: 1"x1" opening

## **ADDITIONAL INFORMATION**

### **Size**

88 x 68 cm x 107- 127 cm (height adjustable)  
33 x 27 inches x 42.5-50 inches