



SKU: S311181

MOLIFT MOVER 180 PATIENT LIFT HOIST

From: \$5,499.99

The Molift Mover 180 Patient Lift Hoist is a lightweight yet durable mobile lifter designed to support everyday transfers in both clinical and home care environments. It offers a 180 kg safe working load and a compact frame for easy manoeuvrability.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Molift Mover 180 Patient Lift Hoist is a reliable, all-purpose patient hoist built for efficient and dignified transfers in aged care, hospital and home environments. Constructed from lightweight aluminium, it balances portability with strength, supporting a safe working load of 180 kg. The ergonomic push bar handle offers multiple grip positions, improving carer comfort and control during operation. Its broad lifting range makes it suitable for various transfers, including floor recovery and high-surface lifts. The Mover 180 comes standard with a four-point sling bar, promoting comfort and optimal positioning.

Molift's eco-friendly NiMH batteries power the hoist, while the built-in service notification system ensures optimal performance. The optional Molift Service Tool allows data extraction for maintenance planning and usage monitoring.

Who it's for

- For clients needing reliable lifting support across floor, bed and seated transfers.
- For therapists wanting a mobile, ergonomic hoist suitable for routine and complex care scenarios.

Key features

- **180 kg safe working load** – supports a wide range of lifting needs
- **Lightweight aluminium frame** – improves manoeuvrability and reduces carer strain
- **Four-point sling bar** – offers stable, comfortable positioning during transfers
- **Ergonomic push bar** – multiple grip positions for easier handling
- **Wide lifting range** – enables lifts from the floor or elevated surfaces
- **Eco-friendly NiMH batteries** – delivers consistent power and reduced environmental impact
- **Service notification system** – alerts carers when servicing is required
- **Molift Service Tool compatible** – enables data tracking for maintenance planning

ADDITIONAL INFORMATION

Brand	Molift
Length	1260 mm
Leg Height	105 mm
Width Chassis, Legs Open Inner	1025 mm
Width Chassis, Outer	670 mm
Castor Size	75 mm
Product Weight	32 kg
Hoisting Speed	50 mm/sec (75 kg)
Hoisting Height Minimum	270 mm
Hoisting Height Maximum	1630 mm
Hoisting Interval	1365 mm
Battery Charge Time	3 Hours
Battery	NiMH 14.4V - 2.6Ah
Rated Performance	50 Hoists at 75 kg
Maximum Weight Capacity	180 kg