




SKU:

ROYALE NORFOLK FITZROY LOWBACK DAY CHAIR




\$799.99

The Royale Norfolk Fitzroy Low Back Day Chair is a robust, adjustable seating solution designed to support individuals with postural or mobility needs. Its bariatric frame and ergonomic features suit a range of care environments.



Scan QR code to
view product online

VARIATIONS

| Image | SKU | Price | Description | Colour |
|--|--------|----------|-------------|-------------------|
|  | S13868 | \$799.99 | | Beige Vinyl |
|  | S13870 | \$799.99 | | Sand Brown Fabric |
|  | S13874 | \$799.99 | | Dark Grey Fabric |

PRODUCT DESCRIPTION

The Royale Norfolk Fitzroy Low Back Day Chair provides a supportive and adjustable solution for those requiring durable bariatric seating. With a contoured low back design and padded seat and armrests, it helps promote posture and comfort during extended periods of sitting.

Height-adjustable legs allow the chair to be tailored to different users, and assist with safer transfers. Its powder-coated steel frame supports a high weight capacity of up to 350 kg.

Available in both vinyl and fabric upholstery, the chair is easy to maintain and adaptable to clinical or home environments.

Who it's for

- For clients who need a low back bariatric chair with personalised height and ergonomic comfort.
- For therapists prescribing a reliable, high-capacity day chair across care settings.

Key features

- **Adjustable seat height** – allows safe and personalised transfer positioning
- **Padded backrest, seat and contoured armrests** – support ergonomic posture and comfort
- **Vinyl or fabric upholstery options** – easy to clean and suited to diverse settings
- **Powder-coated steel frame** – resists corrosion and enhances longevity
- **High weight capacity** – safe working load of up to 350 kg

ADDITIONAL INFORMATION

| | |
|----------------------|----------------|
| Brand | Royale Medical |
| Seat Width | 460mm |
| Seat Depth | 445mm |
| Seat Height | 470-595mm |
| Arm Height | 215mm |
| Back Height | 430mm |
| Overall Width | 585mm |
| Max SWL | 350kg |