




**SKU:** S2314

## **CUP - NOSEY CUP**

**\$16.49**

The nosey cup is ideal for people who have difficulty drinking from standard cups. Features an angled rim so that one can drink without tipping their head. The cup is transparent so a carer can see the level of fluid in relation to the user's mouth. Suitable for people with poor grip, limited neck movement or those who spend long periods in a bed or a wheelchair with a headrest. Has a broad base to reduce tipping, a short stem and tapered body, the cup is easy to grasp and stable. Microwave safe.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

### Key Features

- Microwave and dishwasher safe up to 70°C
- Assists people with poor grip, limited neck movement or those who spend long periods in a bed or a wheelchair with a headrest to drink

## ADDITIONAL INFORMATION

<b>Brand</b>	ILS ADL
<b>Capacity</b>	250 ml
<b>Height</b>	145 mm
<b>Base Diameter</b>	80 mm