




SKU:

OTTOBOCK CLOUD CUSHION

~~\$1,699.99~~ \$1,699.99

The Ottobock Cloud Cushion is a modular seating solution designed for those at high risk of pressure injuries. Its Floam cell construction and customisable design provide exceptional comfort, stability and targeted support, adapting to individual needs.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S16903	\$1,499.99		12" x 12"
	S16898	\$1,499.99		14" x 14"
	S16083	\$1,499.99		14" x 16"
	S7860	\$1,499.99		16" x 16"





Image	SKU	Price	Description	Sizes
	S12720	\$1,469.99		16" x 18"
	S16905	\$1,499.99		17.5" x 16.5"
	S16899	\$1,499.00		17" x 17"
	S16900	\$1,499.00		17" x 18"

Image	SKU	Price	Description	Sizes
	S10710	\$1,499.99		18" x 16"
	S16901	\$1,499.00		18" x 17"
	S9416	\$1,499.99		18" x 18"
	S12724	\$1,699.99		18" x 20"





Image	SKU	Price	Description	Sizes
	S16904	\$1,699.99		19" x 18"
	S12722	\$1,699.99		20" x 16"
	S12721	\$1,699.99		20" x 18"
	S11537	\$1,699.99		20" x 20"





Image	SKU	Price	Description	Sizes
	S10019	\$1,699.99		22" x 18"
	S12723	\$1,699.99		24" x 16"
	S16899	\$1,499.99		Any Sizes
	S16900	\$1,499.99		Any Sizes

Image	SKU	Price	Description	Sizes
	S16901	\$895.99		Any Sizes
	S16902	\$1,699.99		Any Sizes

PRODUCT DESCRIPTION

The Ottobock Cloud Cushion features a modular design that includes up to 30 Floam cells, depending on the cushion size. Floam is a fluid-like polymer that's approximately 30% lighter than traditional fluid materials, delivering effective pressure redistribution and reducing shear. This design allows each Floam cell to be individually configured, supporting postural asymmetries and changing needs. Each cell can be adjusted or replaced as required, ensuring the cushion remains responsive to the individual's needs. A waterproof foam base adds extra protection and durability, while the water-resistant, multi-stretch Dartex outer cover enhances comfort and is easy to maintain.

Who it's for

- For clients at high risk of pressure injuries.
- For therapists seeking a customisable cushion that offers excellent pressure management and postural support for complex seating needs.

Key features

- **Modular Floam cell construction** – allows for individualised pressure redistribution and postural support
- **Lightweight Floam material** – 30% lighter than traditional fluids, reducing overall cushion weight
- **Adjustable cell configuration** – accommodates postural asymmetries and evolving positioning needs
- **Waterproof foam base** – provides additional protection and durability
- **Multi-stretch Dartex cover** – enhances comfort and is water-resistant for easy maintenance

ADDITIONAL INFORMATION

Brand

Otto Bock