




SKU: S822

OVER ARMCHAIR TABLE

\$399.98

The over armchair table has a modern, ergonomically shaped laminated teak table top, which is easy to clean and can be easily tilted to enable users to engage in daily activities such as reading and eating. It comes with a stylish side table to accommodate a mug of coffee/tea which can be fitted to either the left or right hand side of the table.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Over Armchair Table features a modern, ergonomically curved laminated teak table top that is easy to clean and can be tilted for easy use.

ADDITIONAL INFORMATION

Brand	Ausmedic
Max SWL	10 kg (22 lbs)
Main table	109 cm x 41cm 42.9 x 16.1 inches
Side table	20 cm x 20 cm 7.8 x 7.8 inches
Overall width	131 cm (51.5 inches)
Height adjustment	75 cm – 91.5 cm 29.5 – 36 inches
Length between legs	90 cm (35.4 inches)