

SKU: S14681

OVER BED TABLE LARGE

\$299.99

The Royale Medical Over Bed Table is designed for quality and reliability that is unrivalled amongst any other product in its category.

The single lever activated gas strut design makes it easy to adjust and provides versatility and functionality for the user.

The high-quality steel U shape design is durable and fits snugly under the bed to offer a greater quality of life.

The Royale table comes with premium castors which two of which are lockable.”



Scan QR code to
view product online



PRODUCT DESCRIPTION

Key Features:

- Height-adjustable laminated table top. Easy to adjust with a single lever activated gas strut.
- High-quality, long-lasting steel frame
- It comes with four premium castors, two of which are lockable.
- The U-shaped base frame design gives greater stability and allows it to be used with beds and chairs.

ADDITIONAL INFORMATION

Brand	Royale Medical
Total Height	78.5 – 113.5cm
Table Length	79.5cm
Table Width	39.5cm
Base Length	71cm
Base Width	43.5cm
Ground Clearance	10cm
Castors	5cm
Weight	9kg