


SKU: S11199

OVER BED TABLE WITH TILTING SPLIT TABLE

\$399.98

The Royale Medical Over Bed Table with tilting split table is a market leader in design and functionality. Its sleek durable appearance looks stylish in any bedroom setting whilst retaining the versatile functionality of the design.

The table comes in two sections - a flat cornered and stable side compartment with a main adjustable surface perfect for food, reading, drawing, writing and electronic devices at angle. The sturdy height adjustable u-shaped frame is durable and fits snugly under beds with ease. The table comes with four premium castors, two of which are lockable.



Scan QR code to
view product online



PRODUCT DESCRIPTION

Key Features:

- Overbed mobile table with tilting split table top
- Height-adjustable laminated table top
- Easy to adjust with a single lever activated gas strut
- high-quality, long-lasting steel frame
- It comes with four premium castors, two of which are lockable
- U shaped base frame design gives greater stability and allows it to be used with beds and chairs

ADDITIONAL INFORMATION

| | |
|---------------------------|----------------|
| Brand | Royale Medical |
| Side Table Length | 20.5cm |
| Side Table Width | 33.5cm |
| Tilt Table Length | 57cm |
| Tilt Table Width | 38.5cm |
| Total Table Length | 79.5cm |
| Base Length | 80.5cm |
| Base Width | 44.5cm |
| Tilting Slope | 0-68° |
| Max SWL | 125kg |
| Weight | 9.8kg |
| Height | 75cm-109cm |