


**SKU:** S11199

## **OVER BED TABLE WITH TILTING SPLIT TABLE**

**\$399.98**

The Royale Medical Over Bed Table with tilting split table is a market leader in design and functionality. Its sleek durable appearance looks stylish in any bedroom setting whilst retaining the versatile functionality of the design.

The table comes in two sections - a flat cornered and stable side compartment with a main adjustable surface perfect for food, reading, drawing, writing and electronic devices at angle. The sturdy height adjustable u-shaped frame is durable and fits snugly under beds with ease. The table comes with four premium castors, two of which are lockable.



Scan QR code to  
view product online



## **PRODUCT DESCRIPTION**

### Key Features:

- Overbed mobile table with tilting split table top
- Height-adjustable laminated table top
- Easy to adjust with a single lever activated gas strut
- high-quality, long-lasting steel frame
- It comes with four premium castors, two of which are lockable
- U shaped base frame design gives greater stability and allows it to be used with beds and chairs

## ADDITIONAL INFORMATION

<b>Brand</b>	Royale Medical
<b>Side Table Length</b>	20.5cm
<b>Side Table Width</b>	33.5cm
<b>Tilt Table Length</b>	57cm
<b>Tilt Table Width</b>	38.5cm
<b>Total Table Length</b>	79.5cm
<b>Base Length</b>	80.5cm
<b>Base Width</b>	44.5cm
<b>Tilting Slope</b>	0-68°
<b>Max SWL</b>	125kg
<b>Weight</b>	9.8kg
<b>Height</b>	75cm-109cm