


SKU: S11199

OVER BED TABLE WITH TILTING SPLIT TABLE

\$399.98

The Royale Medical Over Bed Table with tilting split table is a market leader in design and functionality. Its sleek durable appearance looks stylish in any bedroom setting whilst retaining the versatile functionality of the design.

The table comes in two sections - a flat cornered and stable side compartment with a main adjustable surface perfect for food, reading, drawing, writing and electronic devices at angle. The sturdy height adjustable u-shaped frame is durable and fits snugly under beds with ease. The table comes with four premium castors, two of which are lockable.



Scan QR code to
view product online



PRODUCT DESCRIPTION

Key Features:

- Overbed mobile table with tilting split table top
- Height-adjustable laminated table top
- Easy to adjust with a single lever activated gas strut
- high-quality, long-lasting steel frame
- It comes with four premium castors, two of which are lockable
- U shaped base frame design gives greater stability and allows it to be used with beds and chairs

ADDITIONAL INFORMATION

Brand	Royale Medical
Side Table Length	20.5cm
Side Table Width	33.5cm
Tilt Table Length	57cm
Tilt Table Width	38.5cm
Total Table Length	79.5cm
Base Length	80.5cm
Base Width	44.5cm
Tilting Slope	0-68°
Max SWL	125kg
Weight	9.8kg
Height	75cm-109cm