


**SKU:** S11199

## **OVER BED TABLE WITH TILTING SPLIT TABLE**

**\$399.98**

The Royale Medical Over Bed Table with tilting split table is a market leader in design and functionality. Its sleek durable appearance looks stylish in any bedroom setting whilst retaining the versatile functionality of the design.

The table comes in two sections - a flat cornered and stable side compartment with a main adjustable surface perfect for food, reading, drawing, writing and electronic devices at angle. The sturdy height adjustable u-shaped frame is durable and fits snugly under beds with ease. The table comes with four premium castors, two of which are lockable.



Scan QR code to  
view product online



## **PRODUCT DESCRIPTION**

### Key Features:

- Overbed mobile table with tilting split table top
- Height-adjustable laminated table top
- Easy to adjust with a single lever activated gas strut
- high-quality, long-lasting steel frame
- It comes with four premium castors, two of which are lockable
- U shaped base frame design gives greater stability and allows it to be used with beds and chairs

## ADDITIONAL INFORMATION

|                           |                |
|---------------------------|----------------|
| <b>Brand</b>              | Royale Medical |
| <b>Side Table Length</b>  | 20.5cm         |
| <b>Side Table Width</b>   | 33.5cm         |
| <b>Tilt Table Length</b>  | 57cm           |
| <b>Tilt Table Width</b>   | 38.5cm         |
| <b>Total Table Length</b> | 79.5cm         |
| <b>Base Length</b>        | 80.5cm         |
| <b>Base Width</b>         | 44.5cm         |
| <b>Tilting Slope</b>      | 0-68°          |
| <b>Max SWL</b>            | 125kg          |
| <b>Weight</b>             | 9.8kg          |
| <b>Height</b>             | 75cm-109cm     |