

SKU: S11193

HERO MEDICAL OVER TOILET AID - ALUMINIUM RUST FREE

\$149.99

The Hero Medical Over Toilet Aid is made from a lightweight aluminium frame. The product has a SWL of 205 kg, height adjustable seat, raised arm supports, and is a durable toileting aid solution. Ideal for users with limited mobility as the raised seat and arm supports can make transferring easier and less strenuous.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Hero over toilet aid has a removable seat and lid for easy cleaning and anti-slip rubber tips to ensure safe bathroom use when surfaces are wet. The product comes with a splash guard and height adjustable extension legs to allow multiple seat height options.

Key Features

- Height adjustable legs
- Removable seat, lid, and splashguard
- Anti-rubber tips
- Raised armrests for additional support and easy transferring on/off
- SWL 205 kg
- Frame has been treated with corrosion resistance and powder coated paint

ADDITIONAL INFORMATION

Brand	Hero Medical
Overall Width	56 cm
Overall Depth	55 cm
Seat Width	34 cm
Width Between the Arms	45 cm
Seat Height Range	41 cm – 54 cm
Armrest Height (from seat to top of arm)	22 cm
Unit Weight	4 kg
Weight Capacity / Max SWL	205 kg