




SKU: S840

OVERBED TABLE - CASTORS

\$169.98

- Slide over beds, chairs, or couches
- Small rails prevent items from sliding off
- Smooth laminate provides a convenient writing surface
- Adjustable in height and comes with swivel castors



Scan QR code to
view product online

PRODUCT DESCRIPTION

- Convenient smooth surface
- It provides a convenient writing, eating, or reading surface and keeps anything you need close by
- The Overbed Table is ideal for those with limited strength and mobility
- Equipped with castors for easy movement
- Features 4 casters that measure 5 cm in diameter
- It also features adjustable height and angle for custom use
- Maximum weight capacity of 10 kg

ADDITIONAL INFORMATION

Brand	Days Healthcare
Max SWL	10 kg
Table size	61 cm x 40 cm
Height	61 cm – 91 cm