



OVERBED TABLE - CASTORS

~~\$165.99~~ \$129.99

Manufactured from chrome plated and powder coated steel with an adjustable laminated teak effect table top. The overbed table is adjustable in height and comes with swivel castors for ease of maneuverability. Built with rails on the edges of the tables to prevent items slipping off.



SKU: S840

Category: [Bedroom](#)

ADDITIONAL INFORMATION

Max SWL	10 kg
Overall Length	610 mm
Overall Width	400 mm
Height Adjustment	635 mm - 910 mm
Brand	Days Healthcare

Manufacturer's part number DAYBED69400