



SKU: S1002

PRIDE PATHRIDER 10 MOBILITY SCOOTER

\$3,149.99

The sleek, sporty Pride Pathrider 10 mobility scooter delivers high performance operation, all-new features and feather-touch disassembly. With a full complement of advanced standard features exclusive low-profile, puncture-proof tyres, back-lit battery gauge, and wraparound easy drive tiller, the Pathrider 10 is a greater value than ever before.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Pride's exclusive one-hand feather-touch dis-assembly provides the ultimate quick transport system. Auto-connecting front-to-rear harness and easy drive tiller with wrap-around handles enhance safety and make steering the Pathrider 10 mobility Scooter easy.

Key Features

- Solid, non-marking tyres
- SWL 159 kg (360lbs)
- Heaviest part is 29.48 kg
- Rugged suspension
- Automatic front and rear lighting system
- Is powered by 2 x 30 AH gel batteries
- Wraparound easy-drive tiller
- Rearview mirror standard
- Front and rear suspension

ADDITIONAL INFORMATION

Brand	Pride
Product Safe Working Load (SWL) - in kg	159 kg
Product Battery Size / Type - in AH / (gel/lead acid/lithium)	2 x 30 AH Lead Acid
Product Range - in km	25 km
Product Speed (Max) - in km/hr	8.5 km/hr
Product Warranty - in months	Mechanical 24 months Electrical 12 months Batteries are covered by a six (6) month warranty.
Product Size - in carton(s) (WxLxH) in cm	58 x 130 x 65
Product Size - in carton(s) (gross weight) in kg	71 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	56.5 x 115 x 88.9
Product Size - overall assembled (gross weight) in kg	93.5 kg
Product Size - seat base width in cm	46 cm
Product Size - seat base depth in cm	43 cm

Product Size - seat back width in cm	45 cm
Product Size - seat back height in cm	35 cm
Product Size - seat to floor in cm	52 cm – 55.8 cm
Product Size - front/ rear tyre type/size	10" / Solid
Product Size - front/rear wheel size in inches	7"
Product Drive - turning circle in cm	157.5 cm
Product Drive - max safe ground clearance in cm	5 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	No
Product Charger Type	3 AH
Product Motor Size	4 Pole / 1344 Watt
Product Lighting	LED
Product Display	Analog