



SKU: S1002

PRIDE PATHRIDER 10 MOBILITY SCOOTER

\$3,149.99

The sleek, sporty Pride Pathrider 10 mobility scooter delivers high performance operation, all-new features and feather-touch disassembly. With a full complement of advanced standard features exclusive low-profile, puncture-proof tyres, back-lit battery gauge, and wraparound easy drive tiller, the Pathrider 10 is a greater value than ever before.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Pride's exclusive one-hand feather-touch dis-assembly provides the ultimate quick transport system. Auto-connecting front-to-rear harness and easy drive tiller with wrap-around handles enhance safety and make steering the Pathrider 10 mobility Scooter easy.

Key Features

- Solid, non-marking tyres
- SWL 159 kg (360lbs)
- Heaviest part is 29.48 kg
- Rugged suspension
- Automatic front and rear lighting system
- Is powered by 2 x 30 AH gel batteries
- Wraparound easy-drive tiller
- Rearview mirror standard
- Front and rear suspension

ADDITIONAL INFORMATION

| | |
|--|--|
| Brand | Pride |
| Product Safe Working Load (SWL) - in kg | 159 kg |
| Product Battery Size / Type - in AH / (gel/lead acid/lithium) | 2 x 30 AH Lead Acid |
| Product Range - in km | 25 km |
| Product Speed (Max) - in km/hr | 8.5 km/hr |
| Product Warranty - in months | Mechanical 24 months Electrical 12 months Batteries are covered by a six (6) month warranty. |
| Product Size - in carton(s) (WxLxH) in cm | 58 x 130 x 65 |
| Product Size - in carton(s) (gross weight) in kg | 71 kg |
| Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm | 56.5 x 115 x 88.9 |
| Product Size - overall assembled (gross weight) in kg | 93.5 kg |
| Product Size - seat base width in cm | 46 cm |
| Product Size - seat base depth in cm | 43 cm |

| | |
|---|--------------------|
| Product Size - seat back width in cm | 45 cm |
| Product Size - seat back height in cm | 35 cm |
| Product Size - seat to floor in cm | 52 cm – 55.8 cm |
| Product Size - front/ rear tyre type/size | 10" / Solid |
| Product Size - front/rear wheel size in inches | 7" |
| Product Drive - turning circle in cm | 157.5 cm |
| Product Drive - max safe ground clearance in cm | 5 cm |
| Product Drive - max safe incline slope in degree | 8° |
| Product Drive - suspension | No |
| Product Charger Type | 3 AH |
| Product Motor Size | 4 Pole / 1344 Watt |
| Product Lighting | LED |
| Product Display | Analog |