


SKU:

INVACARE PEGASUS METRO MOBILITY SCOOTER

\$4,699.99

The Invacare Pegasus Metro focuses on stability and reliability, while equipped with amazing new features to make sure users are safe during each ride. A stable and trusty mobility device, users can enjoy daily trips around town or even take it out for a longer journey. The Pegasus Metro provides independence and efficiency and combines ergonomics with absolute style. You can now enjoy a relaxing drive with the Invacare Pegasus Metro.



Scan QR code to
view product online



VARIATIONS








Image	SKU	Price	Description	Product Colour
	B10028	\$4,699.99		Crystal
	B10024	\$4,699.99		Green
	B10042	\$4,699.99		Chrome
	B10027	\$4,699.99		Red

Image	SKU	Price	Description	Product Colour
	B10025	\$4,699.99		Sapphire Blue
	B10026	\$4,699.99		Silver
	B10023	\$4,699.99		White

PRODUCT DESCRIPTION

The Invacare Pegasus Metro Mobility Scooter is designed with the user's comfort and convenience in mind. It features a comfortable, padded seat that swivels for easy access, and an adjustable steering column that can be positioned to fit the user's preferences. The steering column is also equipped with a digital display that provides real-time information about speed, distance, and battery life.

The scooter also features a range of safety features, including anti-tip wheels, automatic brake system and more. It also has a powerful LED lighting system that provides enhanced visibility in low-light conditions.

Key Features

- **Design:** The Invacare Pegasus Metro scooter has a modern and sleek design, with a stable four-wheeled structure and a quality suspension system that helps to absorb shocks and bumps. It also has a comfortable, padded seat with adjustable armrests, and an adjustable steering column that can be positioned to suit the user's needs.
- **Performance:** The Invacare Pegasus Metro scooter has a maximum speed of 10km/h and a range of up to 54km on a single charge, making it ideal for longer trips and outdoor use. It also has a weight capacity of up to 136kg.
- **Safety features:** The Invacare Pegasus Metro scooter comes with a range of safety features, including anti-tip wheels, a horn, and automatic speed reduction when turning. It also has front and rear lights for improved visibility in low-light conditions.
- **Convenience features:** The Invacare Pegasus Metro scooter includes a number of convenience features, such as a large front basket, a USB charging port, and adjustable seat/armrests and headrest.
- **Shroud Colour options:** The Invacare Pegasus Metro scooter comes in a range of colours.

ADDITIONAL INFORMATION

Model	Pegasus Metro
Weight Capacity / Max SWL	136 kg
Empty Weight	110 kg
Total Length	127 cm
Driving Unit Width	66 cm
Total Height	125 cm
Floor to Seat Height	44 cm – 51 cm
Backrest Height	49 cm
Backrest Width	51 cm
Seat Depth	47 cm
Front/Rear Tyre Type	Pneumatic
Front/Rear Tyre Size	13' inch
Ground Clearance	7.8 cm
Turning Diameter	260 cm
Max Safe Slope	6 degrees
Battery Capacity/Size	2 x 12V/50Ah
Motor Capacity	240w/600w
Charger Type	8AH
Maximum Driving Speed	10Km/h
Maximum Driving Distance/Range	54Km

Carton Size (W x L x H)

80 cm x 153 cm x 75 cm