

SKU:

## INVACARE PEGASUS METRO MOBILITY SCOOTER

**\$4,699.99**

The Invacare Pegasus Metro Mobility Scooter delivers dependable performance with advanced stability, comfort and real-time controls. With a 54 km range and smart speed reduction, it supports independence at every turn.



Scan QR code to  
view product online



**VARIATIONS**





Image	SKU	Price	Description	Product Colour
	B10028	\$4,699.99		Crystal
	B10024	\$4,699.99		Green
	B10042	\$4,699.99		Chrome
	B10027	\$4,699.99		Red

Image	SKU	Price	Description	Product Colour
	B10025	\$4,699.99		Sapphire Blue
	B10026	\$4,699.99		Silver
	B10023	\$4,699.99		White

## PRODUCT DESCRIPTION

The Invacare Pegasus Metro Mobility Scooter is engineered for stable and secure performance across a range of surfaces. It combines long-range capability with a comfort-focused design and a suite of built-in safety features.

A comfortable padded seat offers swivel access, seat depth adjustment, backrest recline and headrest support to suit individual needs.

The adjustable tiller positions easily and includes a digital display that shows real-time speed, distance and battery life.

Safety is a key feature of the Pegasus Metro, with electronic speed reduction around bends, a two-step disengage lever to prevent freewheeling, anti-tip wheels and a full LED lighting system.

Electronics and motor components are also protected from water and corrosion.

Added functionality includes a USB charging port, front carry basket, and traction control for enhanced handling on varied terrain.

The Pegasus Metro has a travel range of up to 54 km.

All Invacare scooters include 12 months of complimentary 24/7 roadside assistance upon registration with the manufacturer.

### Who it's for

- For clients needing a mid-to-large sized scooter with reliable stability and extended range.
- For therapists recommending a smart, safety-focused scooter for independent outdoor use.

### Key features

- **Up to 54 km range with electronic speed reduction** – supports longer trips and safe cornering
- **Padded swivel seat with recline and headrest** – promotes ergonomic comfort
- **LED lighting and brake indicators** – increase visibility in low light
- **Two-step disengage lever and anti-tip wheels** – prevent freewheeling and enhance stability
- **Water and corrosion protected electronics** – ensure long-term durability
- **Traction control button and USB charging port** – adds convenience for daily use
- **Front basket included** – for accessible storage
- **12-month 24/7 roadside assistance** – free with registration through Invacare

## ADDITIONAL INFORMATION

<b>Model</b>	Pegasus Metro
<b>Weight Capacity / Max SWL</b>	136 kg
<b>Empty Weight</b>	110 kg
<b>Total Length</b>	127 cm
<b>Driving Unit Width</b>	66 cm
<b>Total Height</b>	125 cm
<b>Floor to Seat Height</b>	44 cm – 51 cm
<b>Backrest Height</b>	49 cm
<b>Backrest Width</b>	51 cm
<b>Seat Depth</b>	47 cm
<b>Front/Rear Tyre Type</b>	Pneumatic
<b>Front/Rear Tyre Size</b>	13' inch
<b>Ground Clearance</b>	7.8 cm
<b>Turning Diameter</b>	260 cm
<b>Max Safe Slope</b>	6 degrees
<b>Battery Capacity/Size</b>	2 x 12V/50Ah
<b>Motor Capacity</b>	240w/600w
<b>Charger Type</b>	8AH
<b>Maximum Driving Speed</b>	10Km/h
<b>Maximum Driving Distance/Range</b>	54Km

**Carton Size (W x L x H)**

80 cm x 153 cm x 75 cm