




SKU:

PERMOBIL F5 CORPUS VS

The Permobil F5 VS is a standing wheelchair that allows clients to program a standing sequence, providing benefits such as minimising secondary conditions and enabling greater social interaction.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Permobil F5 VS is a highly versatile wheelchair that offers clients the option to stand using a fully programmable standing sequence. With no sacrifices made to functionality or design compared to the regular F5 model, clients can easily transition between sitting and standing positions, which has numerous health benefits. By reducing the risks associated with immobility, such as pressure injuries, contractures, and reduced organ capacity, the F5 VS can significantly improve clients' quality of life. In addition, standing enables clients to perform tasks that would otherwise be impossible and increases social interaction by allowing them to be at the same level as their peers. The F5 VS is an excellent option for anyone seeking a standing wheelchair that offers maximum functionality and customisation.

ADDITIONAL INFORMATION

Brand	Permobil
Drive Wheel Configuration	Front Wheel Drive
Weight Capacity	136 kg
Transit Approved	Yes
Transport Weight	196 kg
Max Speed	10 km/h
Base Width	25.5"
Base Length	43"
Pediatric Seating	No
Seat Width Range	17" - 23"
Seat Depth Range	14" - 22"
Seat To Floor Height Range	17.5" - 19.5"
Posterior Tilt Option	0° - 50°
Anterior Tilt Option	No
Recline Option	85° - 180°
Legrest Elevation Option	85° - 180°
Seat Elevation Option	14"
Power Standing Option	Yes