




SKU:

PERMOBIL K450 MX

The Permobil K450 power chair offers the unique feature of powering the seat to the floor, allowing users to interact with children and access items on the ground. The seating system is adjustable to the user's size and includes a free growth kit.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Permobil K450 power chair is a versatile and unique option that allows users to power the seat to the floor. It is perfect for interacting with children and facilitating transfers. The chair offers a seat-to-floor height range of 3" to 26", allowing users to sit up high at a table or access items on the ground. Additionally, the seating system is completely adjustable and comes with a free growth kit, making it easy to tailor the chair to the size of the user. With these features, the Permobil K450 is a great choice for anyone seeking a customisable and innovative power chair.

ADDITIONAL INFORMATION

Brand	Permobil
Type	Paediatric
Drive Wheel Configuration	Rear Wheel Drive
Weight Capacity	57 kg
Transit Approved	Yes
Transport Weight	152 kg
Max Speed	8 km/h
Base Width	25"
Base Length	35"
Paediatric Seating	Yes
Seat Width Range	10" - 16"
Seat Depth Range	10" - 18"
Seat To Floor Height Range	3" - 26" (Elevator)
Posterior Tilt Option	0° - 45°
Anterior Tilt Option	No
Recline Option	No
Legrest Elevation Option	No
Seat Elevation Option	8"
Power Standing Option	No