


**SKU:** S327151

## **POOL HOIST PORTABLE AQUATIC PAL2 136KG**

**\$17,499.99**



The Pool Hoist Portable Aquatic PAL2 lift is a free-standing pool hoist that is a sleek and modern design with an attractive matte finish. The PAL2 has a 240° rotation to help ensure safe transfers and is easy to move by one person.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

The PAL2 lift is a free standing pool hoist that is a sleek and modern design with an attractive matte finish. It is a portable, free-standing pool lift with a 240° rotation that aids in ensuring a safe transfer.

- 240° rotation helps ensure safe transfer
- Easy to move by one person
- Lift Operator, Intelligent Controller conveniently located on the mast – includes 24 Volt rechargeable battery
- Powder-coated stainless steel and aluminium construction
- Includes: battery, charger, battery console cover, water-resistant hand control, footrest, seatbelt assembly, and armrests
- Lifting capacity: 136 kg

## ADDITIONAL INFORMATION

<b>Shipping Dimensions</b>	102 x 126 x 145 cm (L x W x H)
<b>Weight Capacity</b>	136 kg
<b>Brand</b>	S.R.Smith