




**SKU:** S2098

## **PORTABLE SHOWER TRAY ROYALE**

**\$399.99**

The Shower Tray is a vinyl based unit that makes showering easier. The tray can be placed on top of an existing drain hole or used over an open drain hole, by joining a hose to the side of the passage. The Tray can be used along with a standing shower chair for those with disability problems. The Tray is also easily foldable into a carry bag, allowing quick storage and stress-free portability.

- Comes included with instructions – clearly printed on the carry bag
- Comes included as standard with hose. No need to separately source.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

### Key Features

- Easy to assemble & use
- Carry bag included

## **ADDITIONAL INFORMATION**

<b>Brand</b>	Royale Medical
<b>Dimensions when in use:</b>	Open on floor 1100mm x 1100mm x 75mm Outlet location 250mm x 250mm from corner
<b>Dimensions when stored:</b>	Bag 350mm x 1330mm Bag weight 2kg