




SKU: S2098

PORTABLE SHOWER TRAY ROYALE

\$399.99

The Shower Tray is a vinyl based unit that makes showering easier. The tray can be placed on top of an existing drain hole or used over an open drain hole, by joining a hose to the side of the passage. The Tray can be used along with a standing shower chair for those with disability problems. The Tray is also easily foldable into a carry bag, allowing quick storage and stress-free portability.

- Comes included with instructions - clearly printed on the carry bag
- Comes included as standard with hose. No need to separately source.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Key Features

- Easy to assemble & use
- Carry bag included

ADDITIONAL INFORMATION

Brand	Royale Medical
Dimensions when in use:	Open on floor 1100mm x 1100mm x 75mm Outlet location 250mm x 250mm from corner
Dimensions when stored:	Bag 350mm x 1330mm Bag weight 2kg