


SKU: S7386

PORTABLE WHEELCHAIR RAMP - PERSONAL

\$599.99



Portable, lightweight ramp that allows easier access in and around the home or workplace. Ideal for climbing kerbs or similar obstacles. Passes Australian Standards AS3856-1. Includes carry bag.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Invacare Portable Ramps can be quickly folded into a single compact unit and - allows easy storage or transport in any vehicle or home environment. Created from a unique combination of uni-directional glass and carbon fibre materials.

Ramp Sizes

Length	Width	Weight
70cm	74.5cm	4.7kg

ADDITIONAL INFORMATION

Brand	Invacare
Length	70 cm
Max height capacity	17.5 cm
Ramp width	74.5 cm
Ramp weight	4.7 kg
Safe working load	300 kg