



SKU:

DRIVE MEDICAL PORTLAND RECLINER LIFT CHAIR - TWIN MOTOR


\$1,999.99

The Portland Riser Recliner Lift Chair is designed with comfort in mind. It features an exclusive memory foam seat and its dual motor mechanism allows the user to find their perfect comfort position in minutes.

The Portland Riser Lift Chair Armchair features luxurious Sofia fabric and is available in two colours- Mushroom and Graphite.

Features of the Portland Riser Lift Chair include:

- Reclines to an excellent snooze position.
- High quality upholstery.
- Excellent quiet and smooth actions.
- Built in convenient storage pockets.
- Battery back up facility
- New 2 year warranty
- Weight capacity 150kg



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours	Size
	S11856	\$1,999.99		Graphite Grey	Standard
	S15097	\$1,999.99		Graphite Grey	Petite
	S13069	\$1,999.99		Mushroom	Standard
	S14320	\$1,999.99		Mushroom	Petite

PRODUCT DESCRIPTION

The Portland Riser Recliner Lift Chair is designed with luxurious Sofia Fabric and is equipped with memory foam seat cushions to provide additional comfort. The footrest can be elevated when requires and the backrest can be reclined independently. The lift chair also provides an ergonomically contoured seat and back. To Portland Riser armchair provides a user-friendly four button hand controller.

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Height:	110cm (43.5") ST
	101cm (40") PE
Width:	86cm (34") ST
	83cm (33") PE
Depth:	94cm (37") ST
	90cm (35") PE
Seat Height:	48cm (19") ST
	43cm (17") PE
Seat Depth:	54cm (21.5") ST
	48cm (19") PE
Seat Width:	51cm (20") ST
	47cm (18.5") PE
Max SWL	150KG
Colours	Graphite Grey, Mushroom
Size	Standard, Petite
Brand	Drive Medical