



**SKU:**

## **DRIVE MEDICAL PORTLAND LIFT CHAIR (TWIN MOTOR)**





**\$2,499.99**

The Drive Medical Portland Lift Chair (Twin Motor) provides personalised comfort through independent motor control and a pressure-relieving memory foam seat. Upholstered in soft 'Sofia' fabric, it's available in Petite and Standard sizes.



Scan QR code to  
view product online

**VARIATIONS**

Image	SKU	Price	Description	Colours	Size
	S11856	\$2,499.99		Graphite Grey	Standard
	S15097	\$2,499.99		Graphite Grey	Petite
	S13069	\$2,499.99		Mushroom	Standard
	S14320	\$2,499.99		Mushroom	Petite

## PRODUCT DESCRIPTION

The Drive Medical Portland Lift Chair (Twin Motor) is designed to support individuals who need assistance with sitting and standing transitions. Its twin motor mechanism allows separate adjustment of the backrest and legrest, offering a range of seated and reclined positions for individual comfort and posture needs.

The Portland features a memory foam seat to help distribute pressure evenly, reducing discomfort during extended periods of sitting.

A three-tier waterfall backrest provides enhanced support through the lumbar and upper body. Upholstered in luxurious 'Sofia' fabric in Graphite Grey or Mushroom colours, it combines softness with durability. It is available in both Petite and Standard size options.

For added safety, the hand control includes a built-in key to help prevent unintended operation.

### Who it's for

- For clients who need pressure relief and adjustable recline positioning to aid with standing and sitting.
- For therapists wanting a supportive, memory foam lift chair with independent controls for client comfort.

### Key features

- **Twin motor mechanism** – enables independent control of backrest and legrest
- **Memory foam seat** – helps distribute pressure and increase comfort
- **Three-tier waterfall backrest** – supports spinal alignment and full back comfort
- **Sofia fabric upholstery** – soft, durable and available in two colours
- **Safety key in hand control** – prevents accidental activation
- **Quiet motor operation** – ensures discreet and smooth movement
- **Simple hand control** – easy to operate for daily use
- **Side storage pocket** – keeps essentials within easy reach
- **Snooze-friendly recline** – supports full relaxation
- **Fire retardant materials** – meets required safety standards
- **Available in Petite and Standard sizes** – better fit for different needs

## ADDITIONAL INFORMATION

<b>Height:</b>	110cm (43.5") ST
	101cm (40") PE
<b>Width:</b>	86cm (34") ST
	83cm (33") PE
<b>Depth:</b>	94cm (37") ST
	90cm (35") PE
<b>Seat Height:</b>	48cm (19") ST
	43cm (17") PE
<b>Seat Depth:</b>	54cm (21.5") ST
	48cm (19") PE
<b>Seat Width:</b>	51cm (20") ST
	47cm (18.5") PE
<b>Max SWL</b>	150KG
<b>Colours</b>	Graphite Grey, Mushroom
<b>Size</b>	Petite, Standard
<b>Brand</b>	Drive Medical