



PRIDE MOBILITY POWER GLIDE ASSIST

~~\$1,450.99~~ \$1,100.99

The Power Glide Assist is perfect for medical institutions and caretakers who can easily maneuver wheelchairs with the aid of the Power Glide unit. The Power Glide is designed with a dual wheel option with a maximum speed of 5.5km/h and a range of upto 16km.

The Power Glide Assist features,

- Dual-wheel design for added traction and stability
 - 136 kg (300 lbs.) weight capacity
- Direct drive system provides traction on loose gravel and wet grass
- Maximum speed up to 5.5 km/h (3.4 mph)
 - Range up to 16 km/h (10 miles)
 - Travel bag included
- Extended push handles included
 - Wheelchair is not included

**FREE DELIVERY TO YOUR DOOR.
AUSTRALIA WIDE (ONLINE ONLY)**



SKU: S8484

Categories: [Power Chair Accessories](#),
[Wheelchair Accessories](#)

GALLERY IMAGES



PRODUCT DESCRIPTION

The brand new Power Glide Assist is an excellent unit that aids pushing a manual wheelchair with ease. The Power Glide Assist is perfect for caretakers and medical institutions for daily operational use.

ADDITIONAL INFORMATION

Maximum speed: Up to 5.5 km/h (3.4 mph)

Range per charge:	Up to 16 km (10 miles)
Turning radius:	85 cm (33.8")
Motor:	DC12V, 75W or 280W
Battery:	12V, 18AH
Charger:	1.5 A, off-board
Drive wheel:	20.23 cm (8") solid
Weight capacity:	136 kg (300 lbs.)
Dimensions :	Length: 30 cm (11.81") Width: 20 cm (7.87") Height: 24.5 cm (9.65")
Fits wheelchair widths:	35.5 cm - 50.8 cm (14" - 20")
Total Unit weight:	14 kg (31 lbs.)
Brand	Pride