



SKU:

## SEAT WALKER - 8 INCH PREMIUM BURGUNDY 200KG

~~\$399.99~~ **\$399.99**

The 8-inch Premium Burgundy Seat Walker (200kg) is fitted with 8-inch flat ridged wheels to increase traction and reduce wheel instability. The lightweight aspect to the design and durable structure makes it an ideal indoor and outdoor walker for users.





It also includes a new "light touch" lock brake system with a 200kg safe working limit for the low, medium and large sizes, while the Bariatric size has a 230kg SWL. The frame has been enlarged between the rear wheels to provide more walking space.

The new nylon bag provides 75% more storage space and conceals personal items such as wallets and keys.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Seat
	S11182	\$349.99		High
	S11181	\$349.99		Medium
	S11180	\$349.99		Low
	S11183	\$399.99		Bariatric

## **PRODUCT DESCRIPTION**

This seat walker is built with a heavy-duty frame and wider seat to provide comfort for users up to 200kg (Bariatric 230kg). The 8-inch Premium Burgundy weighs only 8.6kg and is easily collapsed, allowing smooth transporting and manoeuvring.

### **Key Features**

- Heavy duty construction
- Excellent for indoor and outdoor use
- Quality storage space for valuables
- SWL of 200kg for Low, Medium and High seat walkers
- SWL of 230kg for Bariatric model
- Seat height options across the range
- Ergonomic handbrakes for improved safety and stability
- Height adjustable handles
- “Light touch” lock brake system
- Comfortable seat and backrest
- Simple folding mechanism
- Simple installation

## ADDITIONAL INFORMATION

<b>Brand</b>	Royale Medical
<b>Weight Capacity / Max SWL</b>	200kg and 230kg (S11183)
<b>Overall height</b>	Min 76cm Max 87cm Min 87cm Max 99cm Min 91cm Max 102cm
<b>Overall length</b>	70 cm, 70.5 cm, 70.5 cm
<b>Overall width</b>	60 cm
<b>Seat depth</b>	35 cm
<b>Seating width</b>	45 cm Bariatric seat width 55cm
<b>Seat height</b>	50cm, 55cm, 59 cm
<b>Folded width</b>	24.5cm
<b>Castor size</b>	20cm, 8 in
<b>Unit weight</b>	8.6kg – 9kg
<b>Warranty</b>	12 Months Warranty. Please note: Brake maintenance and adjustment is not covered under warranty.