

SKU: S13929

PREMIUM BARIATRIC PADDED SHOWER STOOL - ALUMINIUM RUST FREE

\$699.99

This is the Royale Premium bariatric shower stool with a padded seat. The shower stool has been designed with 'pinched' lightweight aluminium tubing, making it highly durable, functional, and stronger. This product has a safe usage weight of 305 kg and has been tested and complies with Australian Standards. The shower stool has a seating width of 65 cm and has a comfortable padded seat. The Royale premium has height adjustable legs for users to alter the seat height and ensure safety when transferring on and off the seat.



Scan QR code to
view product online



PRODUCT DESCRIPTION

A great showering solution for users needing a wider, heavy-duty seating option. This Royale premium shower stool has been designed with 'pinched' tubing, adding extra strength and stability. This shower chair is heavy-duty and has a maximum usage weight of 305 kg. The stool has a padded seat, anti-slip rubber tips and raised armrests for additional support when lowering and raising. The stool has height adjustable legs which will decrease or increase the height of the seat. All components of the frame have been treated with corrosion resistant powder paint. The product is easily cleaned and maintained.

Key Features

- Lightweight aluminium frame
- Heavy-duty
- 65 cm width between the arms
- Durable anti-slip tips
- Padded seat
- Height adjustable legs
- Raised armrests for stability/support when transferring on and off
- Frame has been treated with corrosion resistant powder paint
- SWL 305 kg

ADDITIONAL INFORMATION

Brand	Royale Medical
Frame	Aluminium
Overall Width	84 cm
Overall Depth	50 cm
Seat Width	50 cm
Width Between the Arms	65 cm
Seat Height Range	43 cm – 55 cm
Armrest Height (from seat to top of arm)	22 cm
Unit Weight	4.8 kg
Weight Capacity / Max SWL	305 kg