


SKU: S11967



PREMIUM BARIATRIC PADDED SHOWER STOOL (56 CM, 280 KG) - ALUMINIUM RUST FREE

\$399.99

This is the Royale Premium bariatric shower stool with a padded seat. The shower stool has been designed with 'pinched' lightweight aluminium tubing, making it highly durable, functional, and stronger. This product has a safe usage weight of 280 kg and has been tested and complies with Australian Standards. The shower stool has a seating width of 56 cm and has a comfortable padded seat. The Royale premium has height adjustable legs for users to alter the seat height and ensure safety when transferring on and off the seat.



Scan QR code to
view product online

PRODUCT DESCRIPTION

A great showering solution for users needing a wider heavy-duty seating option. This Royale premium shower stool has been designed with 'pinched' tubing, adding extra strength and stability. This shower stool is heavy-duty and has a maximum usage weight of 280 kg. The stool has a padded seat, anti-slip rubber tips and raised armrests for additional support when lowering and raising. The stool has height adjustable legs which will decrease or increase the height of the seat. All components of the frame have been treated with corrosion resistant powder paint. The product is easily cleaned and maintained.

Key Features

- Lightweight aluminium frame
- Heavy-duty
- 56 cm width between the arms
- Durable anti-slip tips
- Padded seat
- Height adjustable legs
- Raised armrests for stability/support when transferring on and off
- Frame has been treated with corrosion resistant powder paint
- SWL 280 kg

ADDITIONAL INFORMATION

Brand	Royale Medical
Frame	Aluminium
Overall Width	76 cm
Overall Depth	50 cm
Seat Width	50 cm
Width Between the Arms	56 cm
Seat Height Range	50 cm – 63 cm
Armrest Height (from seat to top of arm)	20 cm
Unit Weight	4.6 kg
Weight Capacity / Max SWL	280 kg