


SKU: S11964



PREMIUM BARIATRIC SHOWER CHAIR - ALUMINIUM RUST FREE

\$399.99

This is the Royale Premium Bariatric Shower Chair. This model is constructed with lightweight 'pinched' aluminium tubing, it is heavy-duty and has a safe working load of 280kg. This shower chair is for users with lower mobility or at risk of falling over whilst showering. This shower chair has the maximum width between the arms for users within our Royale range (56cm) and a high-density plastic seat to minimise flex. Fitted with rubber anti-slip tips to ensure user safety even if surfaces are wet. Compliant with Australian Standards.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Royale Premium Bariatric shower chair is a lightweight showering solution to enable anyone with decreased balance, mobility or strength to shower safely and independently while sitting. The shower chair has height adjustable legs to allow users to adjust the seat height. The raised armrests give users additional support when transferring on and off. One leg extension comes with extra holes which can help stabilise the shower chair if the floor is uneven.

Key Features

- Lightweight aluminium frame
- Bariatric heavy-duty
- Height adjustable legs
- Durable anti-slip tips
- High-density plastic seat with drainage holes
- Frame has been treated with anti-corrosion powder coated paint
- Width between the arms 56cm
- SWL 280kg

ADDITIONAL INFORMATION

Brand	Royale Medical
Frame Material	Aluminium
Overall Width	77cm
Overall Depth	52cm
Seat Width	44cm
Width Between the Arms	56cm
Seat Height Range	44cm – 60cm
Armrest Height (from seat to top of arm)	24cm
Unit Weight	4.1kg
Weight Capacity / Max SWL	280kg