


SKU: S11188

PREMIUM HEAVY DUTY SHOWER CHAIR - ALUMINIUM RUST FREE

\$299.99

The Royale Premium Shower Chair is a heavy-duty option for users. The shower chair is constructed with lightweight 'pinched' aluminium tubing and a high density plastic seat to minimise flex. The chairs 'A' frame provides maximum space between the armrests and offers a safe working load of 280kg. The chair has all the necessary safety features such as durable anti-slip tips, seat drainage holes, raised armrests and Australian standards' certification.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Royale Premium has a heavy-duty 'A' frame shaped design for maximum seat width and comfort. Accompanied by a solid plastic seat with drainage holes, thicker tubing and a durable powder coated finish, makes this shower chair a low maintenance showering solution. Fitted with anti-slip rubber tips, raised armrests and height adjustable legs with a fourth extension leg for perfect balance on uneven surfaces. This is an ideal showering solution to reduce the risk of falls.

Key Features

- Extra wide seating
- Lightweight 'pinched' aluminium tubing
- Heavy duty
- Height adjustable legs
- Durable anti-slip tips
- Raised arm supports for easier transferring on/off
- SWL 280kg
- Australian certified
- Powder coated corrosion resistance paint

ADDITIONAL INFORMATION

Brand	Royale Medical
Frame Material	Aluminium
Overall Width	68cm
Overall Depth	52cm
Seat Width	44cm
Width Between the Arms	50cm
Seat Height Range	49cm – 61cm
Armrest Height (from seat to top of arm)	24cm
Unit Weight	4.7kg
Weight Capacity / Max SWL	280kg