



SKU:

PRESSURE CARE CUSHION - JAY (ION)

~~\$749.99~~ \$749.99

DUAL-LAYERED FOAM BASE WITH
VARYING FIRMNESS

The Ion wheelchair seat cushion features a dual-layered foam base for a combination of structural support and comfort. For optimized immersion for each individual, the foam's firmness is determined by the cushion width, which directly corresponds with client weight.

VISCO™ MEMORY FOAM TOP LAYER

The top layer of soft, Visco memory foam featured in this wheelchair cushion conforms to each individual and distributes weight over a larger area for pressure reduction and maximum comfort.

INNOVATIVE, DUAL-COVER SYSTEM


















For convenience and protection, the JAY Ion features a dual-cover system. The inner cover is water-resistant and easy to wipe off and clean, while its Aqua-guard zipper and anti-wicking thread protect the foam base from fluids. The machine-washable, stretch outer cover reduces surface tension and utilizes naturally anti-microbial, silver-impregnated X-static® fiber. A layer of 3DX spacer fabric promotes airflow and

dissipates heat and moisture for additional skin protection.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S16031	\$499.99		14x15"
	S17397	\$499.99		15x15"
	S16525	\$499.99		15x16"
	S17145	\$499.99		15x18"
		\$449.99		16x16"
	S14928	\$499.99		16x17"
	S15092	\$499.99		16x18"
	S17301	\$499.99		16x19"
	S13270	\$749.99		17x17"
	S13900	\$499.99		17x19"
		\$449.99		18x16"
	S7729	\$499.99		18x17"
	S12237	\$499.99		18x18"
	S11264	\$499.99		18x19"
	S16999	\$499.99		18x20"
	S14397	\$499.99		18x22"
	S14531	\$499.99		20x18"
	S12038	\$499.99		20x20"

ADDITIONAL INFORMATION

Max SWL	136 kg : 14" -21" wide, 225 kg: 22" – 24" wide
Product Weight Max	1.5 kg(based on 18" x 18" size)
Cover Type	X-Static® with Ballistic Nylon
Brand	Sunrise Medical