




SKU:

JAY ION CUSHION

~~\$649.99~~ **\$649.99**

The Jay Ion Cushion provides reliable skin protection and postural stability for people at moderate risk of pressure injury. Its contoured dual-layer foam base and antimicrobial outer cover make it a practical and hygienic solution for daily comfort.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S12670	\$549.99		16" x 16"
	S15092	\$549.99		16" x 18"
	S2410	\$599.99		18" x 16"
	S12237	\$599.99		18" x 18"

Image	SKU	Price	Description	Sizes
	S16999	\$599.99		18" x 20"
	S14397	\$599.99		18" x 22"
	S14531	\$649.99		20" x 18"
	S12038	\$649.99		20" x 20"

PRODUCT DESCRIPTION

The Jay Ion Cushion is designed for active users who need effective skin protection and postural support during daily activities. Its contoured foam base promotes postural alignment and reduces pressure, while the dual-layer construction balances comfort and stability.

A soft Visco memory foam top layer enhances immersion and comfort, while the firmer base provides essential support. The moisture-resistant inner cover features an Aqua-guard zipper and anti-wicking thread to protect the foam from spills and incontinence.

The stretch outer cover includes antimicrobial properties, supporting hygiene in everyday use.

Lightweight and easy to handle, the Jay Ion Cushion suits active individuals needing daily seating support.

Who it's for

- For clients at moderate risk of pressure injury who need a lightweight, comfortable cushion for daily use.
- For therapists prescribing a reliable cushion that balances skin health and postural support for clients in community and home settings.

Key features

- **Contoured dual-layer foam base** – supports posture and reduces pressure during prolonged sitting
- **Soft top Visco memory foam layer** – enhances comfort and immersion
- **Moisture-resistant inner cover** – protects against spills and promotes durability
- **Antimicrobial stretch outer cover** – supports hygiene in daily use
- **Lightweight design** – makes handling and transfers easier for active users

ADDITIONAL INFORMATION

Product Weight Max	1.5 kg(based on 18" x 18" size)
Cover Type	X-Static® with Ballistic Nylon
SWL	100kg to 149kg