

SKU: S20127



## PRIDE PATHRIDER 150XL TURBO MOBILITY SCOOTER WITH 100AH BATTERIES

**From: \$10,299.99**

ILS EXCLUSIVE PATHRIDER 150XL Turbo  
Mobility Scooter

Rugged and built to last, the **Pride Pathrider 150XL Turbo** is everything you want in an outdoor scooter! Equipped with front and rear suspension, this aggressive scooter is designed for cruising on walking trails or driving around your property. Enjoy total visibility with complete LED lighting. With a user-friendly console and dual hydraulic brakes for added safety, the Pathrider 150 XL Turbo is the perfect way to embrace life in the great outdoors!



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

### **YOUR OUTDOOR ADVENTURE AWAITS!**

#### **Key Features**

- Limited recline, black vinyl, high-back seat with lumbar and pillow-top headrest, flip-up height, and angle-adjustable armrests
- CTS Suspension includes adjustable shocks for greater comfort
- Backup sensors with LED light
- Dual hydraulic brakes
- User-friendly LED console displays time, temperature, kilometres driven, and trips taken
- Adjustable delta tiller with ergonomic, wraparound handles / storage compartment included
- USB charger built into the tiller for convenient charging of smart phones and portable devices
- Durable and stylish rear bumper
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Full LED lighting package includes headlights, hazard lights, rear backup sensor LED lights, and directional signals. XLR charger port included
- Tiller shield and basket included

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	159 kg (350 lbs)
<b>Product Battery Size / Type</b>	2 x 100 AH Gel
<b>Product Range</b>	30 km
<b>Product Speed (Max)</b>	10 km/hr
<b>Product Colour</b>	Desert Sand
<b>Product Warranty</b>	Mechanical 24 months Electrical 12 months Batteries are covered by a six (6) month warranty.
<b>Product Size - in carton(s) (WxLxH) in cm</b>	88 x 162 x 87
<b>Product Size - in carton(s) (gross weight) in kg</b>	170 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	83 x 156 x 63.3
<b>Product Size - overall assembled (gross weight) in kg</b>	167.8 kg
<b>Product Size - seat base width in cm</b>	46 cm
<b>Product Size - seat base depth in cm</b>	42 cm

<b>Product Size - seat back height in (inc headrest if applicable) cm</b>	53.5 cm
<b>Product Size - seat to floor (include range if height adjustable) in cm</b>	76 cm – 80 cm
<b>Product Size - front/rear tyre type/size in inches/solid or pneumatic</b>	14.5" / Pneumatic
<b>Product Drive - max safe ground clearance in cm</b>	14 cm
<b>Product Drive - max safe incline slope in degree</b>	10°
<b>Product Drive - suspension</b>	Yes (Front And Rear CTS Suspension)
<b>Product Charger Type</b>	8 AH
<b>Product Motor Size</b>	4 Pole / 2304 Watt
<b>Product Lighting</b>	LED
<b>Product Display</b>	Digital