




SKU:

PRIDE C-101 ELECTRIC RECLINER LIFT CHAIR

\$1,399.99

Elite performance meets understated class. The C-101 power lift recliner from Pride Mobility sets a new standard in customised relaxation. Available in three striking colour options, a durable wooden frame composed of engineered furniture grade laminate/hardwood results in a chair that's strong, supportive, and reliable. This lift chair is single motor and has a SWL of 147 kg.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Suede
	S1871	\$1,399.99		Fern
	S1432	\$1,399.99		Cocoa
	S1872	\$1,399.99		Sandal

PRODUCT DESCRIPTION

The Pride C101 Lift Chair features a comfortable, plush backrest and cushioned seat that are both upholstered in a soft, durable fabric. The chair also includes padded armrests for added comfort and support. With its quiet and smooth motor, this lift chair effortlessly reclines to a comfortable position, and then smoothly lifts to help you stand up from the chair.

- Single motor
- Industry-first hand control with USB charger port
- Quiet & smooth lift system
- Wooden frame composed of engineered furniture grade laminate/hardwood
- Includes a battery backup system
- SWL 147 kg

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Brand	Pride Mobility
Back Width	72 cm
Seat to Floor Height	47 cm
Seat Depth	51 cm
Seat Width	51 cm
Weight Capacity / MAX SWL	147 kg