



SKU:

PRIDE MOBILITY C-101 LIFT CHAIR (SINGLE MOTOR)




\$1,599.99

The Pride Mobility C-101 Lift Chair (Single Motor) supports safe repositioning with single motor functionality. Designed to support individuals up to 147 kg, it combines plush cushioning with intuitive controls.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
	S1871	\$1,599.99		Fern
	S1432	\$1,599.99		Cocoa
	S1872	\$1,599.99		Sandal

PRODUCT DESCRIPTION

The Pride Mobility C-101 Lift Chair (Single Motor) is a medium-sized single motor lift chair that offers everyday support and ease of use.

With its single motor design, the backrest and legrest operate together, creating smooth movement into recline or lift positions.

This model includes a soft fabric finish, padded armrests, and quiet mechanism.

The simple, back-lit hand control features a USB charging port for convenience.

Available in three colour choices – Cocoa, Fern and Sandal.

Who it's for

- For clients who need medium-sized supportive seating with easy lift assistance.
- For therapists wanting a daily-use lift chair with intuitive features for their client.

Key features

- **Single motor operation** – supports simultaneous lift and recline
- **Durable frame and mechanism** – built to last and support up to 147 kg
- **Plush padded armrests and backrest** – for enhanced comfort
- **User-friendly hand control with USB** – for simple control and device charging
- **Battery backup** – ensures function during power interruptions
- **Side pockets** – for easy access to essentials

ADDITIONAL INFORMATION

Back Width	72 cm
Seat to Floor Height	47 cm
Seat Depth	51 cm
Seat Width	51 cm
Weight Capacity / MAX SWL	147 kg