




SKU:

PRIDE MOBILITY C-101 LIFT CHAIR (SINGLE MOTOR)

\$1,599.99

The Pride Mobility C-101 Lift Chair (Single Motor) supports safe repositioning with single motor functionality. Designed to support individuals up to 147 kg, it combines plush cushioning with intuitive controls.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Colours
S1871	\$1,599.99		Fern
S1432	\$1,599.99		Cocoa
S1872	\$1,599.99		Sandal

GALLERY IMAGES



PRODUCT DESCRIPTION

The Pride Mobility C-101 Lift Chair (Single Motor) is a medium-sized single motor lift chair that offers everyday support and ease of use.

With its single motor design, the backrest and legrest operate together, creating smooth movement into recline or lift positions.

This model includes a soft fabric finish, padded armrests, and quiet mechanism.

The simple, back-lit hand control features a USB charging port for convenience.

Available in three colour choices – Cocoa, Fern and Sandal.

Who it's for

- For clients who need medium-sized supportive seating with easy lift assistance.
- For therapists wanting a daily-use lift chair with intuitive features for their client.

Key features

- **Single motor operation** – supports simultaneous lift and recline
- **Durable frame and mechanism** – built to last and support up to 147 kg
 - **Plush padded armrests and backrest** – for enhanced comfort
- **User-friendly hand control with USB** – for simple control and device charging
 - **Battery backup** – ensures function during power interruptions
 - **Side pockets** – for easy access to essentials

ADDITIONAL INFORMATION

Safe Working Load (SWL)	147 kg
Overall Width	82 cm
Overall Depth	75 cm
Overall Height	104 cm
Seat Width	51 cm
Seat Depth	51 cm
Seat to Floor Height	47 cm
Number of Motors	Single
Brand	Pride Mobility

[Pride Mobility Brochure](#)

[Pride Mobility C-101 User Manual](#)

[Pride Mobility Range Comparison Chart](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](#)

[SINGLE VS TWIN VS QUAD MOTOR LIFT CHAIRS - WHAT'S THE DIFFERENCE?](#)