



SKU: S2020

PRIDE C1 ELECTRIC RECLINER LIFT CHAIR - PETITE

\$1,649.99

The Pride C1 Lift Chair is a quiet and smooth three position lift system. An excellent choice for a smaller, affordable motorised riser recliner chair. An integrated emergency battery back-up ensures you stay mobile even in the event of a power shortage.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
	S2020	\$1,649.99		Chocolate
	S301976	\$1,649.99		Merino Beige
	S350	\$1,649.99		Arctic Blue

PRODUCT DESCRIPTION

Enjoy the Pride C1 petite Lift Chair and its great features including a whisper quiet single motor and a steel framework design. The split back design upholstery features fully removable upholstery components and full length, fully padded and upholstered armrests. Available in arctic blue, chocolate or oatmeal fabric upholstery (visit one of our showrooms to view all your options).

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Seat Height	410 mm
Seat Depth	510 mm
Seat Width	460 mm
Overall Width	660 mm
Overall Height	1000 mm
Max SWL	113 kg
Brand	Pride