



SKU:

## PRIDE MOBILITY C6 LIFT CHAIR (TWIN MOTOR)

**\$2,499.99**

The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.



Scan QR code to  
view product online

### VARIATIONS

SKU	Price	Description	Colours
S1874	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Arctic Blue Fabric
S301978	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Merino Beige Fabric
S1878	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Blue Vinyl

SKU	Price	Description	Colours
S0007	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Chocolate Fabric

## PRODUCT DESCRIPTION

The Pride Mobility C6 Lift Chair (Twin Motor) allows separate adjustment of the backrest and legrest, providing greater control for customised positioning.

With a wide 56 cm seat, sewn pillow back and fully removable covers, the chair provides comfort and adaptability.

Available in four colours, the C6 also includes a hand control with USB charging and built-in battery backup for added convenience.

### Who it's for

- For clients who need flexible positioning and custom comfort in a wider chair.
- For therapists wanting independent adjustability and easy-clean design for their clients.

### Key features

- **Twin motor function** – allows separate control of backrest and legrest
  - **Wide 56 cm seat** – for greater comfort and support
  - **Removable upholstery** – supports hygiene and replacement
  - **Durable frame and steel mechanism** – ensures robust performance
- **User-friendly hand control with USB** – easy to operate and charge devices
  - **Battery backup** – maintains use during power outages
  - **Side pockets** – for personal storage

## ADDITIONAL INFORMATION

<b>Safe Working Load (SWL)</b>	136 kg
<b>Overall Width</b>	87 cm
<b>Overall Depth</b>	90 cm
<b>Overall Height</b>	108 cm
<b>Seat Width</b>	56 cm
<b>Seat Depth</b>	53 cm
<b>Seat to Floor Height</b>	48 cm
<b>Number of Motors</b>	Twin
<b>Brand</b>	Pride Mobility

[Pride Mobility Brochure](#)

[Pride Mobility C6 User Manual](#)

[Pride Mobility Range Comparison Chart](#)

## **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

**[SINGLE VS TWIN VS QUAD MOTOR LIFT CHAIRS - WHAT'S THE DIFFERENCE?](#)**