



SKU:

PRIDE MOBILITY C6 LIFT CHAIR (TWIN MOTOR)

\$2,499.99

The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Colours
S1874	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Arctic Blue Fabric
S301978	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Merino Beige Fabric
S1878	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Blue Vinyl

SKU	Price	Description	Colours
S0007	\$2,499.99	The Pride Mobility C6 Lift Chair (Twin Motor) delivers personalised comfort with twin-motor functionality. Users can adjust the backrest and legrest independently for tailored support, while removable upholstery simplifies maintenance.	Chocolate Fabric

PRODUCT DESCRIPTION

The Pride Mobility C6 Lift Chair (Twin Motor) allows separate adjustment of the backrest and legrest, providing greater control for customised positioning.

With a wide 56 cm seat, sewn pillow back and fully removable covers, the chair provides comfort and adaptability.

Available in four colours, the C6 also includes a hand control with USB charging and built-in battery backup for added convenience.

Who it's for

- For clients who need flexible positioning and custom comfort in a wider chair.
- For therapists wanting independent adjustability and easy-clean design for their clients.

Key features

- **Twin motor function** – allows separate control of backrest and legrest
 - **Wide 56 cm seat** – for greater comfort and support
 - **Removable upholstery** – supports hygiene and replacement
 - **Durable frame and steel mechanism** – ensures robust performance
- **User-friendly hand control with USB** – easy to operate and charge devices
 - **Battery backup** – maintains use during power outages
 - **Side pockets** – for personal storage

ADDITIONAL INFORMATION

Safe Working Load (SWL)	136 kg
Overall Width	87 cm
Overall Depth	90 cm
Overall Height	108 cm
Seat Width	56 cm
Seat Depth	53 cm
Seat to Floor Height	48 cm
Number of Motors	Twin
Brand	Pride Mobility

[Pride Mobility Brochure](#)

[Pride Mobility C6 User Manual](#)

[Pride Mobility Range Comparison Chart](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

[SINGLE VS TWIN VS QUAD MOTOR LIFT CHAIRS - WHAT'S THE DIFFERENCE?](#)