

**SKU:** S12886



## **PRIDE GO CHAIR - NEW GENERATION - POWER CHAIR**

**\$3,299.99**

The Pride Go Chair is a compact battery powered wheelchair which disassembles into easy to transport four parts. The Pride Go power chair is rear wheel drive and primarily for indoor use. With solid 200 mm rear wheels and rear anti-tip wheels ensure a safe ride while a moulded plastic seat with padding and synthetic cover with seat adjustable height help to add comfort to your journey.

**[Online Brochure](#)**



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU    | Price      | Description | Colours     |
|----------------------------------------------------------------------------------|--------|------------|-------------|-------------|
|  | S13023 | \$3,299.99 |             | Pearl White |
|  | S12886 | \$3,299.99 |             | Red         |

## **PRODUCT DESCRIPTION**

Pride go chair has a Full length, width adjustable, padded and removable armrests provide additional comfort as needed and a pair of 12 Volt 12 Amp batteries, housed in a box with handle for easy removal enables you to transport the Pride Go Chair, which dismantles in 4 parts.

OMNI-Casters (nylon, spherical-shaped casters) on front and rear to prevent wheel hang-ups

Choice of 6 colors

18" x 17", 5-point swivel seat with latch release

Increased 300 lb. weight capacity

Faster speed - up to 6 kmph

Larger, 18 amp batteries

Improved agility - longer wheel base; wider stance; compact turning radius 25.4"

Dual swing-away, under seat storage bins

Larger, height & width adjustable armrests

Larger footplate

One hand feather-touch disassembly

60" lap belt

## **ADDITIONAL INFORMATION**

|                                |                  |
|--------------------------------|------------------|
| <b>Range</b>                   | Up to 16 km      |
| <b>Speed</b>                   | Up to 5.5 km/h   |
| <b>Backrest Height</b>         | 530 mm           |
| <b>Max SWL</b>                 | 136 kg           |
| <b>Seat Depth</b>              | 430 mm           |
| <b>Seat Width</b>              | 460 mm           |
| <b>Weight of Heaviest Part</b> | 16.5 kg          |
| <b>Brand</b>                   | Pride            |
| <b>Colours</b>                 | Red, Pearl White |