

SKU: S12886




PRIDE GO CHAIR - NEW GENERATION - POWER CHAIR

\$3,299.99



The Pride Go Chair is a compact battery powered wheelchair which disassembles into easy to transport four parts. The Pride Go power chair is rear wheel drive and primarily for indoor use. With solid 200 mm rear wheels and rear anti-tip wheels ensure a safe ride while a moulded plastic seat with padding and synthetic cover with seat adjustable height help to add comfort to your journey.

[Online Brochure](#)



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	S13023	\$3,299.99		Pearl White
	S12886	\$3,299.99		Red

PRODUCT DESCRIPTION

Pride go chair has a Full length, width adjustable, padded and removable armrests provide additional comfort as needed and a pair of 12 Volt 12 Amp batteries, housed in a box with handle for easy removal enables you to transport the Pride Go Chair, which dismantles in 4 parts.

OMNI-Casters (nylon, spherical-shaped casters) on front and rear to prevent wheel hang-ups

Choice of 6 colors

18" x 17", 5-point swivel seat with latch release

Increased 300 lb. weight capacity

Faster speed - up to 6 kmph

Larger, 18 amp batteries

Improved agility - longer wheel base; wider stance; compact turning radius 25.4"

Dual swing-away, under seat storage bins

Larger, height & width adjustable armrests

Larger footplate

One hand feather-touch disassembly

60" lap belt

ADDITIONAL INFORMATION

Range	Up to 16 km
Speed	Up to 5.5 km/h
Backrest Height	530 mm
Max SWL	136 kg
Seat Depth	430 mm
Seat Width	460 mm
Weight of Heaviest Part	16.5 kg
Brand	Pride
Colours	Red, Pearl White