

SKU: S2864

PRIDE MOBILITY GO GO LX WITH CTS SUSPENSION MOBILITY SCOOTER

From: \$2,399.99

The Pride Go Go LX With CTS Suspension combines compact portability with enhanced ride comfort. Featuring specialised suspension, feather-touch disassembly and a 136 kg weight capacity, it supports smooth travel across everyday terrain.



Scan QR code to
view product online

GALLERY IMAGES





PRODUCT DESCRIPTION

The Pride Go Go LX With CTS Suspension is a portable four-wheeled scooter designed for users who want reliable comfort and flexibility. It integrates Comfort-Trac Suspension (CTS) – a leading suspension system that enhances shock absorption over uneven surfaces, delivering a smoother ride in a compact frame.

With feather-touch disassembly, the scooter breaks down into five lightweight pieces for easy transport and storage. The removable battery box offers off-board charging convenience.

The foldable seat has a swivel function and width-adjustable armrests to ensure comfort and easy transfers.

The scooter's ergonomic wraparound delta tiller is adjustable to suit individual needs, while the front-mounted LED headlight illuminates paths in low-light conditions.

Non-marking tyres and a front basket complete this convenient, travel-ready mobility solution.

Who it's for

- For clients needing a compact scooter with enhanced suspension and easy transportability.
- For therapists seeking a lightweight, user-friendly scooter with advanced comfort features.

Key features

- **Feather-touch disassembly into 5 parts** – simplifies transport and storage
- **Comfort-Trac Suspension (CTS)** – improves comfort over uneven terrain
- **Removable battery box with off-board charging** – adds charging flexibility
- **Delta tiller with ergonomic wraparound handles** – supports control and ease of use
- **Foldable seat with swivel and width-adjustable armrests** – promotes access and comfort
 - **LED front headlight** – improves visibility in low light
 - **Non-scuffing tyres and front basket** – ideal for daily errands

ADDITIONAL INFORMATION

Safe Working Load (SWL)	136 kg
Battery Type	Gel
Battery Size	2 x 12V/18 AH
Range	16.9 km
Max Speed	7 km/h
Size (W x L x H assembled)	52 cm x 104 cm x 87cm
Suspension	Front & Rear
Display	Analog
Max Safe Incline Slope	6°
Turning Circle	116.2 cm
Ground Clearance	3.8 cm
Brand	Pride Mobility

[Pride Mobility Go Go LX Brochure](#)

[Pride Mobility Go Go LX User Manual](#)

[Pride Mobility Scooter Specification Chart](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

DOES THE ACID LEAD BATTERY ON THE GOGO LX MOBILITY SCOOTER IMPACT TRAVELLING AIR OR SEA.

GoGO LX Mobility Scooter battery won't impact on air or sea.