

**SKU:** S2864

## **PRIDE GOGO LX MOBILITY SCOOTER WITH CTS SUSPENSION**

**From: \$2,349.99**

The Pride Go-Go LX Mobility Scooter with CTS Suspension (Comfort-Trac Suspension) brings a new level of performance features and value to travel mobility. Featuring CTS independent front and rear suspension, sleek style, feather-touch dis-assembly and standard lighting, the Pride Go-Go LX with CTS Suspension is one of the most advanced travel products available.



Scan QR code to  
view product online



## **PRODUCT DESCRIPTION**

The Go-Go LX 4-wheel travel scooter's advanced design sets a new standard for practicality and comfort. For added stability, the seat post is attached to the front section. Easy Feather Touch™ assembly and disassembly, flat free non scuffing tires, Red and Blue interchangeable color panels, and a front basket are among the standard features you've come to expect from Pride's popular Go-Go lineup.

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	136 kg (299.8 lbs)
<b>Product Battery Size / Type</b>	2 x 17 AH Gel
<b>Product Range</b>	17 km
<b>Product Speed (Max)</b>	6.4 km/hr
<b>Product Colour</b>	Red
<b>Product Warranty</b>	Mechanical 24 months Electrical 12 months Batteries 6 months
<b>Product Size - in carton(s) (WxLxH) in cm</b>	57 x 122 x 65
<b>Product Size - in carton(s) (gross weight) in kg</b>	74 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	52 x 104 x 87
<b>Product Size - overall assembled (gross weight) in kg</b>	53.4 kg
<b>Product Size - seat base width in cm</b>	43 cm
<b>Product Size - seat base depth in cm</b>	43 cm

<b>Product Size - seat back width in cm</b>	43 cm
<b>Product Size - seat back height in cm</b>	34 cm
<b>Product Size - seat to floor in cm</b>	52.7 cm – 56.5 cm
<b>Product Size - front/ rear tyre type/size</b>	7" / Solid
<b>Product Size - front/rear wheel size in inches</b>	4.5"
<b>Product Drive - turning circle in cm</b>	116.2 cm
<b>Product Drive - max safe ground clearance in cm</b>	3.8 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	Yes (Front And Rear CTS Suspension)
<b>Product Charger Type</b>	2 AH
<b>Product Motor Size</b>	2 Pole / 864 Watt
<b>Product Lighting</b>	LED
<b>Product Display</b>	Analog