

SKU: S2864

PRIDE GOGO LX MOBILITY SCOOTER WITH CTS SUSPENSION

From: \$2,349.99

The Pride Go-Go LX Mobility Scooter with CTS Suspension (Comfort-Trac Suspension) brings a new level of performance features and value to travel mobility. Featuring CTS independent front and rear suspension, sleek style, feather-touch dis-assembly and standard lighting, the Pride Go-Go LX with CTS Suspension is one of the most advanced travel products available.



Scan QR code to
view product online



PRODUCT DESCRIPTION

The Go-Go LX 4-wheel travel scooter's advanced design sets a new standard for practicality and comfort. For added stability, the seat post is attached to the front section. Easy Feather Touch™ assembly and disassembly, flat free non scuffing tires, Red and Blue interchangeable color panels, and a front basket are among the standard features you've come to expect from Pride's popular Go-Go lineup.

ADDITIONAL INFORMATION

Product Safe Working Load (SWL)	136 kg (299.8 lbs)
Product Battery Size / Type	2 x 17 AH Gel
Product Range	17 km
Product Speed (Max)	6.4 km/hr
Product Colour	Red
Product Warranty	Mechanical 24 months Electrical 12 months Batteries 6 months
Product Size - in carton(s) (WxLxH) in cm	57 x 122 x 65
Product Size - in carton(s) (gross weight) in kg	74 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	52 x 104 x 87
Product Size - overall assembled (gross weight) in kg	53.4 kg
Product Size - seat base width in cm	43 cm
Product Size - seat base depth in cm	43 cm

Product Size - seat back width in cm	43 cm
Product Size - seat back height in cm	34 cm
Product Size - seat to floor in cm	52.7 cm – 56.5 cm
Product Size - front/ rear tyre type/size	7" / Solid
Product Size - front/rear wheel size in inches	4.5"
Product Drive - turning circle in cm	116.2 cm
Product Drive - max safe ground clearance in cm	3.8 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	Yes (Front And Rear CTS Suspension)
Product Charger Type	2 AH
Product Motor Size	2 Pole / 864 Watt
Product Lighting	LED
Product Display	Analog