



**SKU:** S2174

## **PRIDE JAZZY SELECT 6 POWER CHAIR**

**\$4,499.99**

The Pride Jazzy Select 6 Power Chair features an in-line motor technology for increased efficiency, torque and range coupled with a PG GC 2 Controller. Built with ease of service in mind and easy side access batteries, the Pride Jazzy has been built for ease of use and comfort.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

### **Features**

- 6" castors on front & rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- Articulating rear castor beam
- Easy access charger port and circuit breaker
- Height and depth-adjustable foot platform
- Built with ease of service in mind
- 15.24 cm (6") castors on front and rear for maximum stability

## ADDITIONAL INFORMATION

<b>Brand</b>	Pride
<b>Max SWL</b>	136 kg
<b>Range</b>	Up to 24 km
<b>Speed</b>	6 km/h
<b>Battery Capacity</b>	35 Amp
<b>Seat Depth</b>	500 mm
<b>Max. seating dimensions</b>	46 cm x 51 cm (18" x 20")
<b>Warranty</b>	2 years