



PRIDE JAZZY SELECT 6 POWER CHAIR

~~\$3,600.99~~ \$3,295.99

The Pride Jazzy Select 6 Power Chair features an in-line motor technology for increased efficiency, torque and range coupled with a PG GC 2 Controller. Built with ease of service in mind and easy side access batteries, the Pride Jazzy has been built for ease of use and comfort.

**FREE DELIVERY TO YOUR DOOR
AUSTRALIA WIDE (ONLINE ONLY)**



SKU: S2174

Category: [Power Chairs](#)

ADDITIONAL INFORMATION

Range	Up to 24 km
Battery Capacity	35 Amp

Speed	6 km/h
Max SWL	136 kg
Seat Depth	500 mm
Max. seating dimensions	46 cm x 51 cm (18" x 20")
Brand	Pride
Warranty	2 years